

Table 5. Comparison between single-tube triplex assay, PRV gB assay and virus isolation results using a selected number of clinical samples collected from PRV JS-2012 and PRV Bristol infected pigs.

PRV strain	Pig #	Sample type	PRV triplex assay			PRV gB assay	Virus isolation
			FAM	HEX	Q670	FAM	
JS-2012	17	Brain	25.79	0	23.15	25.9	+
		Tonsil	28.82	0	24.63	29.17	+
		Trigeminal	31.9	0	25.51	31.78	-
	18	Brain	26.99	0	22.3	26.8	+
		Tonsil	22.72	0	25.61	23.4	+
		Trigeminal	28.81	0	28.96	28.72	+
	19	Brain	24.29	0	19.65	23.98	+
		Tonsil	24.35	0	27.27	24.51	+
		Trigeminal	27.32	0	26.99	27.17	+
Bristol	24	Brain	0	23.18	20.77	23.83	+
		Tonsil	0	28.31	26.99	29.48	+
		Trigeminal	0	27.46	28	28.16	+
	25	Brain	0	24.92	20.16	26	+
		Tonsil	0	23.76	25.48	24.86	+
		Trigeminal	0	27.11	26.21	28.01	+
	26	Brain	0	28.94	25.39	29.99	+
		Tonsil	0	20.78	26.87	21.85	+
		Trigeminal	0	30.77	28.14	31.4	-